
A workflow for sequencing, annotating and curating aquatic symbiomes





AQUATIC SYMBIOSIS **GENOMICS**

A workflow for sequencing, annotating and curating aquatic symbiomes

Camilla Eldridge¹, Michael Paulini¹, Caroline Howard¹, Graeme Oathley¹, Radka Platte¹, Catherine McCarthy¹, Nancy Holroyd¹, Nicola Chapman¹, Neha Wali¹, Haoyu Niu¹, Peter Harrison¹, Alexey Sokolov¹, Talal Ibrahim¹, Claudia C. Weber¹, Emmelien Vancaester¹, Marcela Uliano-Silva¹, Ksenia Krasheninnikova¹, Eerik Aunin¹, Shane McCarthy¹, Will Chow¹, Ying Simms¹, James Torrance¹, Alan Tracey¹, Sarah Pelan¹, Joana Collins¹, Damon-Lee Pointon¹, Jonathan Wood¹, Kerstin Howe¹, Michael Sweet², Anne Thompson³, Felipe Porto⁴, Senjie Lin⁴ Patrick Keeling⁵, Helena Villela⁶, Raguel Peixoto⁶, Victoria McKenna¹, Mara Lawniczak¹, Mark Blaxter¹, and members of the ASG consortium (see https://doi.org/10.5281/zenodo.5385621)

The Aquatic Symbiosis Genomics Project

The Aquatic Symbiosis Genomics Project, a collaboration between fourteen international research hubs, is sequencing more than 1000 species from 500 aquatic symbiotic systems to provide insights into complex ecological and evolutionary relationships. The goal is to generate annotated, gold-standard reference genomes for host and symbiont. An initial workflow, from sample preparation to de novo assembly, curation and annotation, has been developed to overcome the challenges associated with large-scale sequencing of diverse aquatic symbiotic systems.

Project status

- 92 species: arrived at Sanger
- 64 species: completed DNA extractions
- 39 species: in sequencing

Sequencing and primary assembly QC data are available on TOLQC https://tolqc.cog.sanger.ac.uk/

Conclusions and future work

- Aquatic symbiome samples have high levels of environmental contamination.
- Ongoing R & D needed for DNA extraction.
- Ongoing development of incorporating CobiontID into the workflow to enable separation of symbionts at read level.




Figure 1, STS and COPO sample management and tracking.



Figure 2. DNA extraction workflow.

Figure 3: The assembly pipeline.



Sample acquisition and management

Samples are collected by partners, snap frozen, and a detailed metadata manifest completed. We sample ethically and legally, including meeting all applicable regulatory compliance steps. These steps are completed before samples are shipped by specialist cold chain couriers in dry ice.


Weighed subsamples are homogenised in lysis buffer/trizol with BioMasher & PowerMasher. A good quality extraction features High Molecular Weight DNA as the main component. The level of RNA degradation is recorded.

Long read data are generated from a low input library protocol followed by sequencing on the PacBio Sequel Ile. For HiC. DNA-protein complexes are cross-linked prior to sample fragmentation and DNA extraction using the Arima-HiC protocol. Libraries are sequenced on Illumina NovaSeq (PE150

Genome assembly and QC

Genomes are assembled from PacBio HiFi reads using HiCanu or Hifiasm and then purged to separate haplotypes. MitoHifi assembles the mitochondrion and eventual plastids. Salsa scaffolds genomes using HiC reads (Arima or Qiagen).

> Fig 4: Read plot of Tetranucleotides coloured by coding density separating mat. from its cobionts



based on taxa using a combination of

pipeline, the CobiontID workflow and

Figure 5: BlobToolKit plot (GC%

an assembly of a fish genome

myxozoan parasite (lower left)

BlobToolkit.

ASG Workflow

Genome Annotation

Curated genomes are annotated using the ASG annotation pipeline that generates a repeat library for each species and uses RNA-seq to train gene model predictors in Braker2. Initial annotations are viewable on the UCSC browser.

Annotated genomes will be hosted by

Research.

Ensembl where final annotations including

protein domain information will be accessible

and a Genome Note describing the assembly will be published in Wellcome Open

Genome curation and postprocessing

Assembled genomes are curated using HiC maps and gEVAL where HiC contact information, coverage and assembly gaps are used in manual curation to scaffold long read assemblies to chromosomal level. Curated genomes are postprocessed and submitted to the ENA. The assembled scaffolds are separated




Figure 8: A dedicated data portal has been setup

providing BioSample display, status tracking and

Figure 7: Viewing annotations




Figure 6. HiC contact map of a colonial sea squirt. Aplidium turbinatum



